

Medical Elizabeth


Introducing Bodyfriend's second medical device.
This product allows you to measure and manage
blood pressure, while receiving a massage.



 **BODYFRIEND**

270° encircle type blood pressure monitor


Conventional blood pressure monitors could measure blood pressure only when the cuff surrounds the wrist a full 360 degrees. Such design harms the usability of the massage chair when applied. Thus, we came up with 270 degree encircle type BP monitor which has never been done before.





Blood pressure treatment massage

Elizabeth Medical is an integrated hardware/software platform that manages health data based on accurate blood pressure measurement. Massage basically activates the parasympathetic nerve, and in addition to this, Elizabeth Medical helps user to control the breathing through the attached speaker with sounds designed to control inhalation/exhalation ratio. This allows the user to directly experience the change in blood pressure value.



Blood pressure management through App

Elizabeth Medical also calculates the likelihood of cardio-cerebrovascular disease such as heart age and stroke based on the measured blood pressure value through the healthcare application systematically combined.

Technical Specifications



59



30


45



Blood Pressure



Blood Pressure Measurement



Link with Mobile App



Brain Massage with Binaural Beat



Advanced XD Massage Module



Thigh Stimulation & Stretching



Advanced Hand & Arm Massage



Magnetic Accupressure for Hands



Full Body Air Massage



Medically Enhanced Massage Modes

Product Name	ELIZABETH Medical	Reclining Angle	115 to 150 degrees
Length	59 inch	Lower Body Angle	0 to 90 degrees
Width	30 inch	Power Supply	AC 110~220V
Height	45 inch	Rated Current	160W
Weight	280 lbs	Power Frequency	50~60Hz

For more information, please visit our VR studio at <https://ces.bodyfriend.com/2021>